



*Model No.: BPM-3DH*

*Model No.: BPM-3FH*



## Microscopes

### Technical Specification:

Head	:	inclind Dual head
Body	:	inclind Binocular Body, Rotatable through 360*
Nose Piece	:	Perfectly centered bearing based, Quadruple revolving nose piece.
Mechanical Stage	:	large stable base (120x 125mm) fitted with co-axial precision graduated mechanical stage for X and Y movement of slide upto 50mm x 75mm
Sub Stage	:	Abbe N.A.1.30 condenser, focussable with rack and pinion through 25mm continously visible iris diaphragm with built in swing out filter holder.
Anti Fungus Objectives	:	Parafoal, Achromatic Din objectives;4x/5x,10x,40x/45x SL, 100x (SL) oil immersion.
Anti Fungus Eye Pieces	:	Highest quality paired Huy.5x / 6x and WF 10x Eye Piece.
Focusing	:	Single Co-axial knob provided based on a four stage gear reduction system for Coarse and fine focusing for silk smooth operation.
Illunimation	:	Fitted with built in base halogen light source ; fitted with 6V-20W Halogen bulbs ; On / Off Switch and light intensity regulator ; complete with power cord and plug. Also a set of plano- concave reflector provided for day light use.

*Same as Model BPM-3FH but fitted with Five Teaching Inclined Head Inclined at 45°*