



## *Student Microscopes*



## **Microscopes** *Student microscopes*

### *FOCUS STUDENT MICROSCOPE*

*MODEL SM-7 ISI MARKED 3686 -1966*

*Triple revolving Nose Piece, Horse Shoe base, Separate Coarse & Fine Motion knobs having stopper for slide safety.*

*ACHROMATIC OBJECTIVE : 10X AND 45X  
HUGENIEN EYE PIECE : H10X AND H15X*

*CONDENSER :*

*Fitted with abbe condenser NA 1.2 and Iris Diaphragm moving up and down with rack and pinion.*

*Packing :*

*Packed in box with dust protection cover with illustrated manual*

### *FOCUS DISSECTING MICROSCOPE*

*MODEL DM-3 ISI MARKED IS: 4228 - 1967*

*EYE PIECE :*

*Aplanatic 10x and Acromatic 20x EP  
on joint arm with glass plate. Mirror for light;  
Rack and Pinion for Focusing*

*PACKING :*

*Packed in box with dust protection cover and  
illustrated manual.*

*Model No.: SM-7*



### *Dissecting Microscope*



*Model No.: DM-3*