



TOOL MAKER
MICROSCOPE



Focus Tool maker's Microscope (T M- 3) : With X-Y Movement on ball bearing slides. Having Micrometer (50mm dia) L.C. 0.01 or 0.005mm With goniometer. Fixed **Magnification 30x** . Having both transmitted and incident light in wo

Optional ACCESSORIES :

- A.Depth Measuring attachment
- B.Centre holding device attachment
- C V - B L O C K
- D V - S U U P O R T
- E.Measuring Eye Pieces with reticles (each)



Abbe Refractometer

Model No.: AR-204



Designed for the Accurate & Quick Examination of the Refractive indices & Mean Dispersion of Liquids, Solids & Powder.

ND Range:

ND 1.3000 to ND 1.7000 in steps of 0.001

Accuracy ± 0.0002

Sugar Percentage.

0 to 95% in steps of 1. Division

Accuracy : $\pm 0.5\%$

Temperature : 0-100°C